


Biopsy – pathology

Tubes with formaldehyde 4%, pH 7.4




- Larger containers upon request
- Container without formaldehyde upon request
- Not suitable for pathogen detection
- Ratio of tissue and formaldehyde: $\geq 1:4$

Unfixed operation preparations, particularly tumor preparations must be transported as quickly as possible to Pathology in order to limit cold ischemia time. Cold ischemia time affects the assessability of certain biological parameters.

Please note the exact sampling time.

Tube 'Michel's Medium'



For unfixed skin biopsies

Shipping material:

Please label with name, first names and date of birth of the patient